

## TROUBLESHOOTING AEROSOL ACCESSORIES

This guide helps verify correct configuration of flexiVent aerosol accessories and troubleshoot issues associated with aerosol generation.

### HARDWARE CONNECTIONS

Make the following connections and ensure that the cables are firmly seated:

- » Insert **Lemo Connector** (red, Figure 1) in **AUX Port 1** (back of flexiVent)
- » Insert **Aeroneb Head Connector** (blue, Figure 1) in **Aeroneb Head** port
- » White cable is well-connected to **Aeroneb Controller**




Figure 1. Aeroneb accessory hardware

**NOTE:** Due to lack of ionic charge, deionized water *does not* nebulize, please use phosphate-buffered saline (PBS).

### AERONEB CONTROLLER STATUS

There are two LED lights on the Aeroneb Controller, please note the following:

- » **Power LED** (red, Figure 2) illuminates with a solid blue light when the flexiVent is powered-on
- » **Nebulization LED** (blue, Figure 2) flashes a blue light *ONLY* when nebulization is active, ie. mist is generated




Figure 2. Aeroneb Controller status LED lights

### SOFTWARE CONSIDERATIONS

Please follow the instructions below:

1. Launch the flexiWare software and enter an Experimentation Session
2. Select an Aeroneb associated template, marked by AN (example *FlexiVent FX Mouse AN* in red, Figure 3)

| Experiment Template Selection  |                   |                   |         |
|--------------------------------|-------------------|-------------------|---------|
| Title                          | Instrument Family | Created           | Status  |
| FlexiVent FX-Mouse Default     | flexiVent FX      | 09/10/20 14:27:24 | Current |
| FlexiVent FX-Mouse OX          | flexiVent FX      | 09/10/20 14:27:29 | Current |
| FlexiVent FX-Mouse AN          | flexiVent FX      | 09/10/20 14:27:34 | Current |
| FlexiVent FX-Mouse AN-OX       | flexiVent FX      | 09/10/20 14:27:40 | Current |
| FlexiVent FX-Mouse AN-EKG      | flexiVent FX      | 09/10/20 14:27:46 | Current |
| FlexiVent FX-Mouse AN-VS       | flexiVent FX      | 09/10/20 14:27:53 | Current |
| FlexiVent FX-Mouse AN-NPFE     | flexiVent FX      | 09/10/20 14:28:01 | Current |
| FlexiVent FX-Mouse AN-NPFE-EKG | flexiVent FX      | 09/10/20 14:28:10 | Current |
| FlexiVent FX-Rat Default       | flexiVent FX      | 09/10/20 14:28:15 | Current |
| FlexiVent FX-Rat AN            | flexiVent FX      | 09/10/20 14:28:21 | Current |
| FlexiVent FX-Rat AN-EKG        | flexiVent FX      | 09/10/20 14:28:27 | Current |
| FlexiVent FX-Rat AN-VS         | flexiVent FX      | 09/10/20 14:28:34 | Current |

Figure 3. Template selection

3. During Aeroneb Priming, hold and visually inspect the Aeroneb Head, note the following:
  - » Aeroneb Head generates aerosol/mist

- » Nebulization LED flashes when mist is generated
4. From an experimentation session, please *right click* in the **Tasks** docker and select *Show All Hidden Tasks* (Figure 4)




Figure 4 . Task Docker – Show All Hidden Tasks

5. Click Nebulizer Task Properties (red, Figure 5)



Figure 5 . Nebulizer – Task Properties

6. In the Task Properties panel, under the *Nebulizer* section, verify that the *Collector Task* is listed as *Dose Response* (see Figure 6).




Figure 6 . Collector Task – Dose Response

7. In the Task Properties panel, under the *Digital Outputs* section, verify that *Aeroneb* is listed and checked, below *Active Digital Outputs* (see Figure 7)




Figure 7. Digital Output – Aeroneb

**NOTE:** Should nebulization issues persist, please contact support for additional recommendations.

## CONTACT INFORMATION

### **North America**

TechSupport@scireq.com

Tel: 1 (514) 286-1429 option 2

Toll Free: 1 (877) 572-4737 option 2

### **Europe & Middle East**

sales@emka.fr

Tel: +33 (1) 40 60 76 00

### **China**

info@bjgyd.com

Tel: +86 (0)10-85376382

### **Japan**

info@emkatech.jp

Tel: +81 (0)6-6476-7135